
sartorius stedim
biotech

Directions for Use | Gebrauchsanleitung | Instruction d'utilisation
Instrucciones para el uso | Istruzioni per l'uso

Sartolab[®]-P20 | Sartolab[®]-P20 plus

83052-000-50

Contents

1. Product Description and Applications
2. Ordering Information
3. Notes
4. Specifications for Sartolab P20
5. Specifications for Sartolab P20 Plus
6. Directions for Use
7. Chemical Compatibility

1. Product Description and Applications

The Sartolab P20 and Sartolab P20 Plus ready-to-use filtration units, for single use only, are delivered sterile and individually packaged.

They have been developed for fast, biocompatible sterile filtration of aqueous solutions and media containing serum.

A surfactant-free cellulose acetate membrane with very low non-specific protein binding characteristics is sealed in a polycarbonate housing.

The combination of the hydrophobic PTFE membrane for automatic venting and the hydrophilic cellulose acetate membrane results in a very high flow rate during filtration.

The glass fiber prefilter (which is 100% free of binding agents) incorporated in the Sartolab P20 Plus may even double the total filtration volume.

2. Ordering Information

The units are available in various configurations:
with either a hose nipple (barb) or luer lock as the inlet connector.
The order numbers for a package of 10 units are listed below:

Sartolab® P20 and Sartolab® P20 Plus

	Sartolab P20		Sartolab P20 plus	
Order number	180 52D	180 53D	180 56D	180 58D
Cellulose acetate membrane 0.2 µm	+	+	+	+
Glass fiber prefilter			+	+
Filling bell	+	+	+	+
Hose nipple inlet	+		+	
Luer Lock inlet		+		+
Hose nipple outlet	+	+	+	+

3. Notes

- The Sartolab P units are constructed so that filtration may only be performed in one direction. The inlet side is labeled “IN” and has eight vent openings arranged around the inner circle. The outlet side is labeled “OUT.”
- When handling the units, make sure to prevent contamination of the filtrate outlet, or the mounted filling bell.
- If the filtrate flow rate drops below the required rate (the filter is blocked), replace the Sartolab P with a new one.
- Sartolab P units are for single use only – do not re-use them.

4. Specifications for the Sartolab P20

Order numbers	180 52D, 180 53D
Housing material	Polycarbonate
Filter material	Surfactant-free cellulose acetate 0.2 µm pore size
Filter diameter	64 mm
Filtration area	20 cm ²
Hold-up volume before bubble point	1 ml
Connector: inlet	Female luer lock or stepped hose nipple with 6–12 mm outer diameter
Connector: outlet	Hose nipple (see above)
Flow direction	Unidirectional
Housing burst pressure	> 5 bar (72 psi)
Bubble point	≥ 3.2 bar (46 psi)*
Max. recommended inlet pressure	3 bar (43 psi)
Flow rate for water	300 ml/min at Δp = 1 bar (14.5 psi)
Filtration range	100 ml – max. 5 l

5. Specifications for the Sartolab P20 Plus

Order numbers	180 56D, 180 58D
Housing material	Polycarbonate
Filter material	Surfactant-free cellulose acetate, 0.2 µm pore size Glass fiber, 100% free of binding agents
Filter diameter	64 mm
Filtration area	20 cm ²
Hold-up volume before bubble point	Approx. 1.5 ml
Connector: inlet	Female luer lock or stepped hose nipple with 6–12 mm outlet diameter
Connector: outlet	Hose nipple (see above)
Flow direction	Unidirectional
Housing burst pressure	> 5 bar (72 psi)
Bubble point	≥ 3.2 bar (46 psi)*
Max. recommended inlet pressure	3 bar (43 psi)
Flow rate for water	300 ml/min at Δp = 1 bar (14.5 psi)
Filtration range	100 ml – max. 10 l

* ≥ 3.1 bar when measured with an automatic integrity tester

6. Directions for Use

6.1. The Sartolab P20 and Sartolab P20 Plus are pressure filtration units.

The liquid to be filtered can be supplied via a pressure tank or a peristaltic pump. Especially for this range of applications, we offer a self-regulating diaphragm pump with electronic speed control and a set of connecting hoses for the Sartolab P20 units (18059). Please read the operating instructions provided with the pump.

6.2. After you have set up the pressure source and filtrate receptacle, remove a Sartolab P20 unit from the box and make sure that the packaging is not damaged.

6.3. Make sure that the inlet side of the pressure filtration unit is located below the seam for opening the package (marked with a “Y”). The inlet side can be identified by the eight vent openings arranged around the inner circle (see Figure 1). The sides of the unit are also labeled “IN” = inlet and “OUT” = outlet. If you have ordered Sartolab P20/Sartolab P20 Plus with a filling bell, it is on the outlet side.

Once you have properly positioned the unit, tear open the packaging at the designated seam (see Figure 2).

6.4. Connect the hose from your pressure source to the inlet side of the Sartolab P20/Sartolab P20 Plus (luer lock connector or stepped hose nipple) without taking the sterile side out of the packaging. Check that the hose is connected tightly enough to withstand the pressure applied.

6.5. Without touching the sterile side, remove the filtration unit from the packaging and mount the unit above the filtrate receptacle, for instance, using a stand (see Figure 3).

6.6. Remove the lid from the filling bell before beginning filtration.

6.7. Start filtration by slowly increasing the applied pressure. The entrapped air escapes through the PTFE membrane by automatic venting. As soon as the cellulose acetate membrane has completely wetted out, you can increase the inlet pressure up to 3 bar.

6.8. After the filtration run has been completed, you can check the integrity of the filter by connecting a pressure source to the unit. Adjust the inlet pressure to a constant 3.0 bar (300 kPa or 44 psi) and hold the Sartolab filtration unit under water in a beaker.

If the filter is intact, no air bubbles should escape from the outlet (downstream) side. Air bubbles that escape through the automatic vents on the upper part do not have any effect on the integrity of the surfactant-free cellulose acetate membrane.

6.9. After completing filtration, discard the filtration unit. If applicable, follow the current regulations for radioactive or infectious material.

Caution

When filtering solutions containing surface-active substances (organic solvent components, alcohols), the water penetration pressure of the PTFE membrane may decrease so much that liquid seeps out of the automatic vents. If necessary, rinse these substances out of the PTFE membrane first, then dry the unit in a drying oven overnight at 105°C.

Avoid unnecessary air intake during the filtration by ensuring, for example, that the inlet end of the intake tubing is kept below the surface of the liquid being filtered. Otherwise, particularly when using Sartolab P20 Plus, a pocket of air may build up above the membrane, resulting in partial or complete “blocking” of the membrane. Should this occur, however, pump all of your solution out of the system first, then continue the filtration run as described in step 6.7.

7. Chemical Compatibility

As mentioned above, the Sartolab P20/ Sartolab P20 Plus is designed for filtration of aqueous solutions – not for organic solvents, alcohols, acetone, surfactants (emulgators), strong acids or bases. For such purposes, we recommend the Sartorius Stedim Biotech Midisart 2000 units, containing a PTFE membrane, order no. 17804E (0.45 µm pore size) or 17805E (0.2 µm pore size).

In the interest of the further development of Sartorius Stedim Biotech products, we reserve the right to make changes to the design or specifications without notice.

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Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part of this
publication may be reprinted or
translated in any form or by any
means without the prior written
permission of Sartorius Stedim
Biotech GmbH.

The status of the information,
specifications and illustrations
in this manual is indicated by
the date given below. Sartorius
Stedim Biotech GmbH reserves
the right to make changes to
the technology, features, specifi-
cations, and design of the
equipment without notice.

Status:
March 2008,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany