

FV-2400 Micro-Spin, Mini-Centrifuge/Vortex




DESCRIPTION

Mini-centrifuge/vortex Micro-Spin **FV-2400** is specially designed for genetic engineering research (for PCR-diagnostics experiments). Units can be used in microbiological, biochemical, clinical laboratories and industrial biotechnological laboratories. Micro-Spin provides simultaneous mixing and separation of samples, using centrifuge and mixing modules, located on the common spin-module.

FV-2400 is an "open type" centrifuge (without lid), that increases the speed of centrifugation and resuspension operations.

SPECIFICATION

Continuous and impulse operation modes	+
Overall dimensions (WxDxH)	120x170x120 mm
Weight	1.4 kg
Nominal operating voltage	120 or 230 V; 50 Hz
Rotation speed (fixed)	2800 rpm
Max. RCF	500 x g
Power consumption (230V / 120 V)	30 W (0.13 A) / 30 W (0.27 A)
----	----
Nominal operating voltage	120 or 230 V; 60 Hz
Rotation speed (fixed)	3500 rpm
Max. RCF	700 x g
Power consumption (230V / 120 V)	25 W (0.1 A) / 30 W (0.27 A)



CAT. NR.

	FV-2400 incl. rotors R-1.5M, R-0.5/0.2M & power plug
BS-010201-AAA	230VAC 50/60Hz Euro plug
BS-010201-AAB	230VAC 50/60Hz UK plug
BS-010201-AA3	230VAC 50/60Hz AU plug
BS-010201-AAC	100VAC 50/60Hz US plug
BS-010201-AAC	120VAC 60Hz US plug
-----	Optional accessories
BS-010205-DK	Rotor R-2/0.5/0.2
BS-010205-CK	Rotor R-2/0.5
BS-010202-AK	Rotor SR-16
BS-010201-EK	Rotor SR-64

FV-2400 Micro-Spin, Mini-Centrifuge/Vortex



BS-010201-CK	IQ OQ document
BS-010201-DK	PQ document
-----	Replacement parts
BS-010201-AK	Rotor R-1.5M
BS-010201-BK	Rotor R-0.5/0.2M

ACCESSORIES




R-1.5M
rotor for 12 x 1.5/2 ml microtest tubes



R-0.5/0.2M
rotor for 12 x 0.5 ml and 12 x 0.2 ml microtest tubes



R-2/0.5
rotor for 8 x 2/1.5 ml and 8 x 0.5 ml microtest tubes




R-2/0.5/0.2
rotor for 6x2/1.5 ml + 6x0.5 ml + 6x0.2 ml microtest tubes

FV-2400 Micro-Spin, Mini-Centrifuge/Vortex



SR-16

rotor for 2 x 8-section 0,2 ml microtube strips



SR-64

rotor for 8 x 8-section 0,2 ml microtube strips - for any type of strips including paired