

Safety Data Sheet

SECTION 1: Identification of the substance/mixture and of the company/undertaking

Identification of the substance or mixture

Product code 10328011
Product name RNASE AWAY

Company/undertaking identification

Life Technologies Corporation
5781 Van Allen Way
PO Box 6482
Carlsbad, CA 92008
+1 760 603 7200

Life Technologies
5250 Mainway Drive
Burlington, ONT
CANADA L7L 6A4
800/263-6236

24 hour Emergency Response for Hazardous Materials [or Dangerous Goods] Incident. Spill, Leak, Fire, Exposure, or Accident. Call CHEMTREC Within the USA + Canada: 1-800-424-9300 and 1-703-527-3887
Outside the USA + Canada: 1-703-741-5970

Country Specific Emergency Number (if available):

For research use only. Not for use in diagnostic procedures

SECTION 2: Hazards identification

GHS Classification

Signal Word
WARNING

Hazard pictograms



Health hazards

Revision date 14-Apr-2022
Product code 10328011

Page 1 / 9
Product name RNASE AWAY

Skin corrosion/irritation	Category 2
Serious eye damage/eye irritation	Category 2

Physical hazards

Not Hazardous

Environmental hazards

Not Hazardous

Hazard Statements

H315 - Causes skin irritation

H319 - Causes serious eye irritation

Precautionary Statements

Prevention

P280 - Wear protective gloves/protective clothing/eye protection/face protection

P264 - Wash hands thoroughly after handling

Response

P302 + P352 - IF ON SKIN: Wash with plenty of soap and water

P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

P332 + P313 - If skin irritation occurs: Get medical advice/attention

P337 + P313 - If eye irritation persists: Get medical advice/attention

P362 + P364 - Take off contaminated clothing and wash it before reuse

Storage

Not Applicable

Disposal

Not Applicable

Other hazards

Not Applicable

HMIS

Health	2
Flammability	0
Reactivity	0

SECTION 3: Composition/information on ingredients

Chemical Name	CAS No.	Common name	EINECS-No	Weight-%
Alkali hydroxide	proprietary	Not applicable	-	<0.1

We recommend handling all chemicals with caution.

SECTION 4: First aid measures

Description of first aid measures

Skin contact	Wash off immediately with plenty of water for at least 15 minutes. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Immediate medical attention is required.
Eye contact	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Immediate medical attention is required.
Ingestion	Never give anything by mouth to an unconscious person. Do not induce vomiting without medical advice. Get medical attention if symptoms occur.
Inhalation	Remove to fresh air. If not breathing, give artificial respiration. If symptoms persist, call a physician.
Notes to Physician	Treat symptomatically.

Most important symptoms and effects, both acute and delayed

H315 - Causes skin irritation H319 - Causes serious eye irritation

Indication of any immediate medical attention and special treatment needed

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SECTION 5: Firefighting measures

Extinguishing media

Suitable extinguishing media	Water spray. Carbon dioxide (CO ₂). Foam. Dry chemical.
Unsuitable extinguishing media	No information available.

Special hazards arising from the substance or mixture

No data

Advice for firefighters

Wear self-contained breathing apparatus and protective suit.

SECTION 6: Accidental release measures

Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation
Avoid contact with skin, eyes or clothing
Use personal protection equipment
See section 8 for more information

Environmental precautions

No special environmental precautions required. Avoid discharge into drains and waterways whenever possible.

Methods and material for containment and cleaning up

Soak up with inert absorbent material.

Reference to other sections

See section 8 for more information.

SECTION 7: Handling and storage

Precautions for safe handling

Always wear recommended Personal Protective Equipment. See section 8 for more information. Do not get in eyes, on skin, or on clothing. Do not ingest. If during normal use the material presents a respiratory hazard, use only with adequate ventilation or wear appropriate respirator.

Conditions for safe storage, including any incompatibilities

Keep in a dry, cool and well-ventilated place. Keep in properly labeled containers. Store in accordance with local regulations.

Specific end use(s)

For research use only.

SECTION 8: Exposure controls/personal protection

Control parameters

Chemical Name	OSHA PEL	OSHA PEL (Ceiling)	ACGIH OEL (TWA)	ACGIH OEL (STEL)
Alkali hydroxide	None	None	None	None

Chemical Name	Brazil - OEL - TWAs (LTs)	Brazil - OEL - Ceilings	Brazil - OEL - Skin Designations
Alkali hydroxide	None	None	None

Engineering Measures Ensure adequate ventilation, especially in confined areas.

Exposure controls

Personal Protective Equipment

Respiratory protection In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

Hand protection Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.

Eye protection Tight sealing safety goggles.

Skin and Body Protection Wear laboratory coat for body protection.

Hygiene Measures Handle in accordance with good industrial hygiene and safety practice

Environmental exposure controls

No special environmental precautions required. Avoid discharge into drains and waterways whenever possible.

SECTION 9: Physical and chemical properties

Information on basic physical and chemical properties

Physical state	liquid	
Color	No data	
Odor	No data	
Molecular Weight	No data	
Melting point / melting range	°C No data	°F No data
Boiling point / boiling range	°C No data	°F No data
Flammability (solid, gas)	No data	
Lower explosion limit	No data	
Upper explosion limit	No data	
Flash point	°C No data	°F No data
Autoignition Temperature	°C No data	°F No data
Decomposition temperature	°C No data	°F No data
pH	10-11.5	
Evaporation rate	No data	
Viscosity	No data	
Solubility	No data	
Partition coefficient: n-octanol/water	No data	
Vapor Pressure	No data	
Specific gravity	No data	
Relative density	No data	
Vapor density	No data	
Explosive properties	No data	
Oxidizing properties	No data	
Particle characteristics	No data	

Other information

Information with regard to physical hazard classes

No information available

Other safety characteristics

No information available

SECTION 10: Stability and reactivity

Reactivity	None known.
Chemical stability	Stable under normal conditions.
Possibility of hazardous reactions	Hazardous reaction has not been reported.
Conditions to avoid	None known.
Incompatible materials	No dangerous reaction known under conditions of normal use.
Hazardous decomposition products	No known hazardous decomposition products.

SECTION 11: Toxicological information

Information on toxicological effects

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Alkali hydroxide	No data available	No data available	No data available

Principal Routes of Exposure

Acute toxicity	Data are conclusive but insufficient for classification.
Skin corrosion/irritation	Skin irritation
Serious eye damage/irritation	Causes serious eye damage
Respiratory or skin sensitization	Data are conclusive but insufficient for classification
Specific target organ toxicity (STOT) – single exposure	Data are conclusive but insufficient for classification
Specific target organ toxicity (STOT) – repeated exposure	Data are conclusive but insufficient for classification
Carcinogenicity	Data are conclusive but insufficient for classification
Germ cell mutagenicity	Data are conclusive but insufficient for classification
Reproductive toxicity	Data are conclusive but insufficient for classification
Aspiration hazard	Data are conclusive but insufficient for classification

SECTION 12: Ecological information

Toxicity

The environmental impact of this product has not been fully investigated.

Chemical Name	Toxicity to algae	Toxicity to daphnia and other aquatic invertebrates	Toxicity to fish	Microtox Data	log Pow
Alkali hydroxide	No data available	No data available	No data available	No data available	No data available

Mobility in soil No information available.

Persistence and degradability No information available.

Bioaccumulative potential No information available.

Results of PBT and vPvB assessment

No information available.

Other adverse effects No information available.

SECTION 13: Disposal considerations

Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations.

SECTION 14: Transport information

IATA / ADR / DOT-US / IMDG

Not regulated in the meaning of transport regulations

UN number or ID number	Not Applicable
UN proper shipping name	Not Applicable
Transport hazard class(es)	Not Applicable
Packing group	Not Applicable

Marine pollutant

No

Special precautions for user

Not Applicable

Maritime transport in bulk according to IMO instruments

Not Applicable.

SECTION 15: Regulatory information

Component	US TSCA
Alkali hydroxide proprietary (<0.1)	-

US Federal Regulations

SARA 313

This product is not regulated by SARA.

Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product does not contain HAPs

US State Regulations

Revision date 14-Apr-2022
Product code 10328011

Page 7 / 9
Product name RNASE AWAY

California Proposition 65

This product does not contain any Proposition 65 chemicals.

WHMIS Hazard Class

D2B - Toxic materials

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the SDS contains all the information required by the CPR.

National Regulations - Brazil

Chemical Name	CAS No.	Brazil - National Agency for Sanitary Surveillance (ANVISA)	Brazil - National List of Carcinogen Agents to Humans (LINACH)
Alkali hydroxide	proprietary	Not Listed	Not Listed

SECTION 16: Other information

Reason for revision SDS sections updated.
Revision number 5
Revision date 14-Apr-2022

For research use only. Not for use in diagnostic procedures.

References

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

Abbreviations and acronyms

TWA - Time-Weighted Average
OELs - Occupational Exposure Limits
STEL - Short Term Exposure Limit
DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List
KECL - Korean Existing and Evaluated Chemical Substances
ENCS - Japan Existing and New Chemical Substances
IECSC - China Inventory of Existing Chemical Substances
PICCS - Philippines Inventory of Chemicals and Chemical Substances
AICS - Australian Inventory of Chemical Substances
NZIoC - New Zealand Inventory of Chemicals
EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances
TSCA - United States Toxic Substances Control Act Section 8(b) Inventory
CEPA - Canadian Environmental Protection Act
EPA - Environmental Protection Agency
OSHA - Occupational Safety and Health Administration of the US Department of Labor
IATA - International Air Transport Association
DOT - Department of Transportation
IMDG - International Maritime Dangerous Goods
ACGIH - American Conference of Governmental Industrial Hygienists
NIOSH - National Institute for Occupational Safety and Health
AIHA - American Industrial Hygiene Association
HMIS - Department of Defense Hazardous Materials Information System
NTP - National Toxicology Program
IARC - International Agency for Research on Cancer

Revision date 14-Apr-2022
Product code 10328011

Page 8 / 9
Product name RNASE AWAY

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"

End of Safety Data Sheet