

41 TD4 Low Speed Automatic-balance Centrifuge



Features:

- ✧ Microprocessor control, DC brushless motor offers stable and low noise operation.
- ✧ Parameters can be set freely during running.
- ✧ Time counts down in seconds when the centrifugation time is less than 1min.
- ✧ Self-diagnosis system can automatically detect the faults of over-speed, no speed signal, and the door lid interlock system during running, ensuring the operation safety.
- ✧ Compact design and small footprint offers ideal performance.

Technical Parameters:

Model	TD4
Order No.	0205194000
Max. Speed	4000rpm
Max. RCF	1980xg
Max. Capacity	8 x 15ml
Speed Accuracy	±30rpm
Time Setting Range	1min to 99min
Noise	<62dB(A)
Power Supply	AC220V±22V 50/60Hz 2A

Total Power	100W
Dimensions (W x D x H)	320x370x235mm
Package Size (W x D x H)	530x410x290mm
Net Weight	11kg
Rotor Optional:	
	Order No: 0301809030 Max. Speed : 4000rpm Max. RCF: 1980xg Max. Capacity: 8 x 15ml