

Conical Tube

(#50015/51015/50050/50515/50250/50040/50115/51115/50150/51150)

1. Description

Conical Tubes are widely used in laboratories from storage of samples to high speed centrifugation. SPL Life Sciences provides very high strength and non-toxic conical tubes to meet the demands.

2. Features

- ❖ Tight sealing screw cap
- ❖ External graduations with marking area
- ❖ Packed in zipper bags
- ❖ Polystyrene (transparent) conical tubes (Cat. No. 51015, 51115, 51150)
- ❖ Skirted, self-standing type (Cat. No. 50515, 50250)
- ❖ For high speed centrifugation (Cat. No. 50040)
- ❖ Provided in racks (Cat. No. 50115, 51115, 50150, 51150)

3. Functional Test

Sterility	:	Passed
Non-cytotoxic	:	Passed
Non-pyrogenic	:	Passed
DNase/RNase - free	:	Passed
Human DNA - free	:	Passed

4. Product Range

In Sleeves

Cat.No.	Bottom Type	Material (Body/Cap)	External Dimension (mm)	Total Volume (ml)	RCF Rating	Sterile	Packaging
50015	Conical	PP / HDPE	17.00 x 120.00	15.00	13,000 x g	+	50 in sleeve / 500
51015	Conical	PS / HDPE	17.00 x 120.00	15.00	3,000 x g	+	50 in sleeve / 500
50050	Conical	PP / HDPE	30.00 x 115.00	50.00	14,000 x g	+	25 in sleeve / 500
50515	Skirted	PP / HDPE	17.00 x 120.00	15.00	-	+	50 in sleeve / 500
50250	Skirted	PP / HDPE	30.00 x 116.00	50.00	-	+	25 in sleeve / 500
50040	Conical	PP / HDPE	30.00 x 115.00	40.00	20,000 x g	+	10 in sleeve / 100

In Rack

Cat.No.	Bottom Type	Material (Body/Cap)	External Dimension (mm)	Total Volume (ml)	RCF Rating	Sterile	Packaging
50115	Conical	PP / HDPE	17.00 x 120.00	15.00	13,000 x g	+	25 in rack / 500
51115	Conical	PS / HDPE	17.00 x 120.00	15.00	3,000 x g	+	25 in rack / 500
50150	Conical	PP / HDPE	30.00 x 115.00	50.00	14,000 x g	+	25 in rack / 300
51150	Conical	PS / HDPE	30.00 x 115.00	50.00	3,000 x g	+	25 in rack / 300

SPL-TDS-COTU

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 4, November 24th, 2023

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.


6. Specifications & Technical Drawing

(1) 15 ml Conical Tube (# 50015/51015/50515/50115/51115)

- Specifications

	Dimensions	Material (Resin)	Packaging
Cap	Diameter: 21.30 mm Height: 11.00 mm	High-density polyethylene (HDPE)	<ul style="list-style-type: none"> ▪ Pieces/sleeve: 25 or 50 ▪ Pieces/case: 500
Body	Diameter: 16.90 mm Height: 119.20 mm Volume: 15.00 ml	Polypropylene(PP) or Polystyrene(PS)	

- Technical Drawing




(2) 50 ml Conical Tube (# 50050/50040/50250/50150/51150)

- Specifications

	Dimensions	Material (Resin)	Packaging
Cap	Diameter: 34.00 mm Height: 13.00 mm	High-density polyethylene (HDPE)	<ul style="list-style-type: none"> ▪ Pieces/sleeve: 10 or 25 ▪ Pieces/case: 100 or 300 or 500
Body	Diameter: 29.30 mm Height: 113.90 mm Volume: 50.00 ml	Polypropylene(PP) or Polystyrene(PS)	

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-COTU

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 4, November 24th, 2023