

Stacker Microtube Rack

(#61048)

1. Description

SPL Stacker Microtube Racks are uniquely designed for handling and storage of Microcentrifuge Tubes.

2. Features

- ❖ Stackable feature for space saving
- ❖ 5 colors: Pink, Yellow, Green, Orange, Blue
- ❖ Non-autoclavable


3. Product Range

Cat. No.	Material	Description	Packaging
61048	ABS	6 x 8 (48 Holes)	10

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs




6. Specifications & Technical Drawing

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 123.00 mm Width: 183.00 mm Height: 26.00 mm	ABS	▪ Pieces/Case: 10

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-STTURA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 20th, 2016