

Cryo Box

(#80025/80081/80181/80281)

1. Description

Cryo Boxes are useful for storing vials in ultra low temperature environments without distortion of storage boxes.

2. Features

- ❖ For storing of 1.5 ml microtubes or cryogenic vials
- ❖ Designed for enhanced drainage of liquid nitrogen
- ❖ Alphanumeric indications for sample location
- ❖ Temperature range: -196 ~ +121°C (Cat. No. 80025, 80081, 80281)
- ❖ 3 colors: Blue, Pink, Green (Cat. No. 80025, 80081, 80281)
- ❖ Cardboard Box: Moisture repellent coating (Cat. No. 80181)
- ❖ Cat. No. 80025, 80081, 80181: 1.2 ml / 1.8 ml Cryovials are available
- ❖ Cat. No. 80281: 5 ml Cryovial use only

3. Product Range

Cat. No.	Material	Feature	Description	Packaging
80025	PC	3 Colors	5 x 5 (25 holes)	1 / 20
80081	PC	3 Colors	9 x 9 (81 holes)	1 / 10
80281	PC	3 Colors	9 x 9 (81 holes)	1 / 6
80181	Cardboard	-	9 x 9 (81 holes)	40

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs



SPL-TDS-CRB

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 3, November 20th, 2018



6. Specifications & Technical Drawing


(1) Cryo Box – PC (# 80025, 80081, 80281) – Specifications Dimensions

– Specifications


		Dimensions		Material (Resin)
Cat. No. 80025		Cat. No. 80081		
Length:	76.83 mm	Length:	133.50 mm	Polycarbonate (PC)
Width:	76.83 mm	Width:	133.50 mm	
Height:	35.17 mm	Height:	35.20 mm	
		Length:	133.50 mm	
		Width:	133.50 mm	
		Height:	96.00 mm	

Packaging		
Cat. No. 80025	Cat. No. 80081	Cat. No. 80281
<ul style="list-style-type: none"> ▪ Pieces/Sleeve: 1 ▪ Pieces/Case: 20 	<ul style="list-style-type: none"> ▪ Pieces/Sleeve: 1 ▪ Pieces/Case: 10 	<ul style="list-style-type: none"> ▪ Pieces/Sleeve: 1 ▪ Pieces/Case: 6

- Technical Drawing



Cat. No. 80025



Cat. No. 80081 / 80281

SPL-TDS-CRB

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 3, November 20th, 2018

(2) Cryo Box – (# 80181) – Specifications Dimensions

– Specifications

Dimensions		Material (Resin)	Packaging
Cap	Body		
Length:	132.5 mm	Cardboard	▪ Pieces/Case: 40
Width:	132.5 mm		
Height:	35.0 mm		
Length:	127.4 mm		
Width:	127.4 mm		
Height:	50.0 mm		

www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr