

Storage Box

(#80100/81100/81081)

1. Description

SPL Storage Boxes are designed to hold microtubes and vials in wide temperature ranges. Two different sizes are available for various purposes.

2. Features

- ❖ Temp. range: -70 ~ +140°C
- ❖ 5 colors: Blue, Green, Orange, Pink, Purple (Cat. No. 81100)
- ❖ Hinged lid type (Cat. No. 80100, 81100)
- ❖ Autoclavable

3. Product Range

Cat. No.	Material	Color	Description	Packaging
80100	PP	Natural Color	10 x 10 (100 Holes)	30
81100	PP	5 colors	10 x 10 (100 Holes)	30
81081	PP	Natural Color	9 x 9 (81 Holes)	40

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs



SPL-TDS-STBO

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 20th, 2016


6. Specifications & Technical Drawing

- Specifications


Dimensions		Material (Resin)	Packaging	
Cat. No. 81081	Cat. No. 80100/81100		Cat. No. 81081	Cat. No. 80100/81100
Length: 130.50 mm	Length: 140.50 mm	Polypropylene (PP)	▪ Pieces/Case: 40	▪ Pieces/Case: 30
Width: 130.50 mm	Width: 140.50 mm			
Height: 35.20 mm	Height: 30.00 mm			

- Technical Drawing

Cat. No. 81081



Cat. No. 80100/81100



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-STBO

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 20th, 2016