

MOUZIP®

(# 82101)

1. Description

MOUZIP® is a unique device that serves as a complete testing bench for a variety of mouse experiments, composed of an experimental mat absorbing all liquids that may occur during experiments and a paper disposal case that wraps the mouse cadaver for ethical handling, storing and disposing.

- The top side of the mat efficiently absorbs any excess liquids during the experiment, while the waterproof bottom maintains hygienic conditions.
- The paper case connected to the mat is folded up in a triangular coffin-like structure to enable ethical storage and hygienic disposal of cadaver after experiments.

2. Features

- ❖ The top side of the mat efficiently absorbs any excess liquids during the experiment, while the waterproof bottom maintains hygienic conditions.
- ❖ The paper case connected to the mat is folded up in a triangular coffin-like structure to enable ethical storage and hygienic disposal of cadaver after experiments.

3. Product Range

Cat.No.	Material (Outer / Inner)	Cover Dimension (cm)	Mat Dimension (cm)	Packaging
82101	Kraft / Fiber mat	13.00 x 24.00	25.00 x 30.00	50

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photograph


6. Specifications & Technical Drawing

(1) MOUZIP® (# 82101)

- Specifications

Dimensions		Material (Outer / Inner)	Packaging
Outer	Inner		
Width - 13.00 cm Length - 24.00 cm	Width - 25.00 cm Length - 30.00 cm	Kraft / Fiber mat	50 per case

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
 To place an order, contact your local distributor or
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-82101

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
 © SPL Life Sciences Co., Ltd., Version 3, July 21st, 2016