

Loop & Needle (#90001/90010)

1. Description

Disposable Loop and Needle consists of a loop on one side and a needle on the other, used for the inoculation of microorganisms. Two different sizes are available.

2. Features

- ❖ For microbe inoculation
- ❖ Dual purpose : a loop on one end and needle on the other
- ❖ Color coded for sizes (White for 1 $\mu\ell$, Blue for 10 $\mu\ell$)

3. Functional Test

Sterilized : Passed
RNase/DNase free : Passed
DNA free : Passed

4. Product Range

Cat. No.	Material	Size ($\mu\ell$)	Length (mm)	Sterile	Packaging
90001	PP	1.00	195.00	+	10 / 500
90010	PP	10.00	195.00	+	10 / 500

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

6. Product Photographs



SPL-TDS-LONE

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 3, November 22th, 2018


7. Specifications & Technical Drawing

- Specifications


Dimensions		Material (Resin)	Packaging	
			90001	90010
Length mm)	195.00	Polypropylene (PP)	<ul style="list-style-type: none"> ▪ Pieces/Sleeve: 10 ▪ Pieces/Case: 500 	<ul style="list-style-type: none"> ▪ Pieces/Sleeve: 10 ▪ Pieces/Case: 500

- Technical Drawing

Cat. No. 90001



Cat. No. 90010



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
 To place an order, contact your local distributor or
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-LONE

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
 © SPL Life Sciences Co., Ltd., Version 3, November 22th, 2018