

CHAPTER 17

Centrifuge Tubes


All of our sterile centrifuge tubes are manufactured out of the highest grade resins and are certified to be RNase, DNase, DNA, pyrogen, and ATP free.

Polypropylene tubes are autoclavable to 121°C at 15 P.S.I. for 15 minutes, and are freezable to -80°C.

- 15,000g centrifuge force for 50mL & 15mL.
- 25,000g centrifuge force for micro centrifuge tube.
- Leak-proof seals.
- Racked & Bulk pack available, recyclable plastic rack reduces waste in the lab.
- Sterile is an option.




50mL Centrifuge Tube, PP, conical bottom, screw cap, 15000g, sterile



Part No.	Volume	Package Type	Pack
CT50B	50mL	bulk pack	25/bag, 20bags/case
CT50R	50mL	rack pack	25/rack, 20racks/case

15mL Centrifuge Tube, PP, conical bottom, screw cap, 15000g, sterile




Part No.	Volume	Package Type	Pack
CT15B	15ml	bulk pack	25/bag, 20bags/case
CT15R	15ml	rack pack	25/rack, 20racks/case

Micrcentrifuge Tube, pp, clear, graduated, 25000g



Part No.	Volume	Package Type	Pack
MCT50	5.0ml	bulk pack	200/bag, 10bags/case
MCT20	2.0ml	bulk pack	500/bag, 10bags/case
MCT15	1.5ml	bulk pack	500/bag, 10bags/case
MCT05	0.5ml	bulk pack	1000/bag, 10bags/case

50mL Centrifuge Tube, PP, conical bottom with skirt, self-standing, screw cap, 3500g, sterile



Part No.	Volume	Package Type	Pack
CT50RS	50ml	rack pack	25/rack, 20racks/case
CT50BS	50ml	bulk pack	25/bag, 20bags/case