


4

CHAPTER


Culture Dish

- Four specifications are available: 35 mm, 60 mm, 100 mm and 150 mm.
- Adherent cell culture and suspended cell culture are available.
- The raised outer edge of the dish cover ensures the stable stacking during transportation.
- Excellent transparency and flatness, more suitable for microscope observation.
- The bottom of the dish is marked to facilitate the localization of cells under the microscope.
- The hand-held ring design makes it easy to wear gloves and grasp.




35mm Cell Culture Dish, PS, Sterile

| Part No. | Volume | Culture area | TC | Specification | Pack |
|----------|--------|--------------------|----|---------------|---------------------|
| CL-D035 | 3ml | 8.8cm ² | ✖ | 35*10mm | 20/bag, 25bags/case |
| CL-DT035 | 3ml | 8.8cm ² | ✓ | 35*10mm | 20/bag, 25bags/case |




60mm Cell Culture Dish, PS, Sterile

| Part No. | Volume | Culture area | TC | Specification | Pack |
|----------|--------|---------------------|----|---------------|---------------------|
| CL-D060 | 5ml | 21.5cm ² | ✖ | 60*15mm | 10/bag, 50bags/case |
| CL-DT060 | 5ml | 21.5cm ² | ✓ | 60*15mm | 10/bag, 50bags/case |



100mm Cell Culture Dish, PS, Sterile

| Part No. | Volume | Culture area | TC | Specification | Pack |
|----------|--------|---------------------|----|---------------|---------------------|
| CL-D100 | 12.5ml | 56.7cm ² | ✖ | 100*20mm | 10/bag, 30bags/case |
| CL-DT100 | 12.5ml | 56.7cm ² | ✓ | 100*20mm | 10/bag, 30bags/case |



150mm Cell Culture Dish, PS, Sterile

| Part No. | Volume | Culture area | TC | Specification | Pack |
|----------|--------|--------------------|----|---------------|--------------------|
| CL-D150 | 35ml | 145cm ² | ✖ | 150*20mm | 5/bag, 20bags/case |
| CL-DT150 | 35ml | 145cm ² | ✓ | 150*20mm | 5/bag, 20bags/case |



Petri Dish, PS

| Part No. | Diam | Pack |
|----------|------|----------------------|
| PD90 | 90mm | 20/bag, 25bags/case |
| PD60 | 60mm | 20/bag, 50bags/case |
| PD35 | 35mm | 10/bag, 230bags/case |